

# Current sources modules for non-invasive modulation of electrophysiology

Alexander Hunold, Ethan Rich, Sebastian Berkes, Max Breuer, Patrick Stein, Robert Guggenberger, Klaus Schellhorn  
neuroConn GmbH, Ilmenau, Germany

## Takeaways

Here we introduce **current source modules** optimally designed for **low-noise transcranial** electrical stimulation or **high-power, high-bandwidth** electrical stimulation with **precise temporal control** and **synchronization** enabling combined **transcranial and neuromuscular stimulation**.

## Background

Investigations of **direct effects of electrical stimulation** on electrophysiological activity requires **low noise stimulation** enabling **synchronous EEG** [1] or **in vitro cell activity** [2] recordings.

Application scenarios of electrical stimulation **target area located deeper in the brain** [3] and include **peripheral nerve** and **muscle stimulation** [4], which require **higher intensities**.

High **temporal precision** is required for both, **temporal interference** stimulation and **paired transcranial and neuromuscular stimulation**.

Stimulation effects on electrophysiological activity need to be assessed in **real-time**

## Innovation

**Current source modules** designed for

**transcranial** stimulation based on **low-noise operational amplifiers**

**peripheral** stimulation based on **high-power, high-bandwidth transistor** gain stages


**modular design** with

**flexible control** via pulse width modulated (**PWM**) signal,

**galvanic insulation** via digital isolators,

**temporal synchronization** via **real-time bus**

## Low-noise current source



3 db edge frequency: **5 kHz**

**100 kHz**

full scale rise time: **1 μs**

**Dynamic range**

**+/- 4 mA** arbitrary waveform

up to 1 ms duration: **+/- 40 mA**






low frequency waves: **+/- 15 mA**

**Key performance**

**Low 1/f noise** for EEG combination

**fast high intensity** for periphery

## Wideband universal isolated power current source



## Conclusion

**new current source modules** enable:

**central** and **peripheral** electrophysiological modulation with

extended parameter sets including

**amplitude modulation** via PWM control

to address further targets



**temporal synchronization** allows:

**pairing of transcranial** or **peripheral** stimulation for

investigations of corticospinal **interactions** and realization of

**temporal interference** stimulation

## Temporal synchronization and galvanic isolation enable temporal interference stimulation



## References

- [1] Wunder, S., Hunold, A., Fiedler, P., Schlegelmilch, F., Schellhorn, K., & Haueisen, J. (2018). Novel bifunctional cap for simultaneous electroencephalography and transcranial electrical stimulation. *Scientific reports*, 8(1), 1-11..
- [2] Ahtiainen, A., Leydolph, L., Tanskanen, J., Hunold, A., Haueisen, J., Hyttinen, J. (2023). Establishing and assessing temporal interference electrical stimulation in neuronal cultures in vitro. *SfN Annual Meeting*.
- [3] Shan, Y., Wang, H., Yang, Y., Wang, J., Zhao, W., Huang, Y., ... & Zhao, G. (2023). Evidence of a large current of transcranial alternating current stimulation directly to deep brain regions. *Molecular Psychiatry*, 1-9.
- [4] Guggenberger, R., Schmidberger, L., Schill, L., & Gharabaghi, A. (2023). Phase-dependent modulation of human corticospinal plasticity by associative pairing of transcranial and neuromuscular stimulation. *bioRxiv*, 2023-04..