

Effectiveness of Transcranial magnetic stimulation

rTMS



We are here for you!

neurocare group

Albert-Einstein-Straße 3, 80331 Ilmenau, Germany

phone: +49 (0)3677 68979-0

fax: +49 (0)3677 68979-15

e-mail: info@neurocaregroup.com

web: www.neurocaregroup.com

rTMS Overview 05

Depression 06

Obsessive compulsive disorder (OCD) 09

Anxiety 10

Post-traumatic stress disorder (PTSD) 11

Pain 12

Stroke 14

Aphasia 16

Multiple sclerosis 17



Overview

Evidence-based guidelines on the therapeutic use of repetitive transcranial magnetic stimulation (rTMS): An update (2014–2018)

AUTHORS	Lefaucheur JP, Aleman A, Baeken C, Benninger DH, Brunelin J, Di Lazzaro V, Filipovic SR, Grefkes C, Hasan A, Hummel FC, Jääskeläinen SK, Langguth B, Leocani L, Londero A, Nardone R, Nguyen JP, Nyffeler T, Oliveira-Maia AJ, Oliviero A, Padberg F, Palm U, Paulus W, Poulet E, Quartarone A, Rachid F, Rektorová I, Rossi S, Sahlsten H, Schecklmann M, Szekely D, Ziemann U
SOURCE	Clinical Neurophysiology, volume 131 (2020), pp 474–528
OBJECTIVE	Overview of the latest research results on the therapeutic use of rTMS and evaluation of the scientific evidence.
METHODOLOGY	Classification of studies in terms of criteria randomized/controlled, number of patients, clear statements on procedure, primary outcome, exclusion criteria, statistical analysis, patient groups. Summary of outcomes.

RESULTS

Level A	Depression (high-frequency rTMS) Neuropathic pain Post-stroke contralateral hand motor recovery in post-acute stage
Level B	Depression (low-frequency rTMS, bilateral rTMS) Post-traumatic stress disorder Fibromyalgia Parkinson’s disease: depressive symptoms Parkinson’s disease: motor symptoms Post-stroke ipsilateral hand motor recovery in post-acute stage Post-stroke nonfluent aphasia in chronic stage Lower limb spasticity in multiple sclerosis
Level C	Obsessive compulsive disorder Addiction/craving Schizophrenia: auditory hallucinations & negative symptoms Tinnitus Hemispatial neglect in post-acute stage after stroke Epilepsy Alzheimer’s disease: cognitive function Complex Regional Pain Syndrome - type I

CONCLUSIONS	A review of hundreds of scientific studies has shown that rTMS can provide significant improvement in various neurological and psychiatric disorders. rTMS is scientifically proven to be definitely effective in the treatment of neuropathic pain and depression as well as in the hand motor recovery from stroke. Furthermore, scientific studies show that rTMS is probably effective for patients with motor impairment, Parkinson’s disease, multiple sclerosis, post-traumatic stress disorder or aphasia.
--------------------	--


Depression

Simultaneous rTMS and psychotherapy in major depressive disorder: Clinical outcomes and predictors from a large naturalistic study

AUTHORS	Donse L, Padberg F, Sack AT, Rush AJ, Arns M
SOURCE	Brain Stimulation 11 (2018) 337-345
OBJECTIVE	Evaluation of feasibility and clinical outcome of rTMS in combination with psychotherapy in major depressive disorder. Identification of clinical predictors of response and remission.
METHODOLOGY	Naturalistic study including 196 patients with depression, rTMS treatment with simultaneous psychotherapy. 97% had been unsuccessfully treated with medication at least once before.

PROTOCOL

Target area	left DLPFC / right DLPFC
Stim. frequency	10 Hz / 1 Hz
Stim. intensity	110 - 120 % RMT
Number of pulses per session	1,500 / 1,200
Treatment	2 - 10 per week
Number of sessions	21 on average
Assessment	BDI, DASS, response defined as $\geq 50\%$ reduction in BDI score



RESULTS

Assessment	
Response	66.3%
Remission	56.0 %
BDI	baseline 31.3, last session 14.1, 55,9 % reduction, $p < 0.001$
DASS D	baseline 28.6, last session 12.1
DASS A	baseline 13.7, last session 5.6
DASS S	baseline 22.4, last session 10.1
Follow-up 6 months Responder group	
BDI	baseline 29.4, last session 8.0, follow-up 13.8
DASS Depression	baseline 27.0, last session 6.5, follow-up 11.8
DASS Anxiety	baseline 14.1, last session 4.0, follow-up 7.2
DASS Stress	baseline 22.3, last session 7.3, follow-up 12.5

CONCLUSIONS

The combined therapy rTMS + psychotherapy led to high response and remission rates. In addition, the effect also showed good stability in the follow-up. Considering the high proportion of treatment resistance against antidepressants in the sample, the result is of clinical relevance.

Depression

Efficacy and Safety of Transcranial Magnetic Stimulation in the Acute Treatment of Major Depression: A Multisite Randomized Controlled Trial

AUTHORS O'Reardon JP, Solvason HB, Janicak PG, Sampson S, Isenberg KE, Nahas Z, McDonald WM, Avery D, Fitzgerald PB, Loo C, Demitrack MA, George MS, Sackheim HA


SOURCE Biol Psychiatry 62 (2007) 1208–1216

OBJECTIVE Effectiveness and safety of TMS over left DLPFC in the acute treatment of major depression.

METHODOLOGY Double-blind, placebo-controlled study with 301 patients with depression (155 active rTMS, 146 sham stimulation). Patients were medication-free and had shown medication resistance prior to rTMS treatment.

PROTOCOL

Target area	Left DLPFC using the 5 cm rule
Stim. frequency	10 Hz
Stim. intensity	120 % RMT
Number of pulses per session	3,000
Treatment	daily, 5 times per week
Number of sessions	max. 30
Assessment	MADRS, HAMD17, HAMD24



RESULTS

Assessment	Significance active vs. sham	Active group	Sham group
MADRS response rate*	P < 0.01	23.9 %	12.3 %
HAMD24 response rate	P < 0.05	23.9 %	15.1 %
MADRS remission rate**	P < 0.05	14.2 %	5.5 %
HAMD24 remission rate***	P < 0.05	17.4 %	8.2 %
Adverse events			
Scalp discomfort/pain reported as mild to moderate		35.8 % diminished after first 5 sessions	3.8 %
Increase in suicidality		1 event	10 events

*50 % improvement from baseline **total score < 10 ***total score < 8 ****total score < 11

CONCLUSIONS

rTMS administered over the left DLPFC using the parameters reported here for a period of up to 6 weeks was effective in treating major depression with a good tolerability profile. These results indicate that rTMS offers clinicians a novel alternative in the treatment of this disorder.

Depression

Efficacy of prefrontal theta-burst stimulation in refractory depression: a randomized sham-controlled study.

AUTHORS Li CT, Chen MH, Juan CH, Huang HH, Chen LF, Hsieh JC, Tu PC, Bai YM, Tsai SJ, Lee YC, Su TP


SOURCE Brain 137 (2014) 2088–2098

OBJECTIVE Antidepressant effect of right pre-frontal continuous, left-prefrontal intermittent and combined theta-burst stimulation compared to sham group of patients showing different levels of antidepressant treatment resistance.

METHODOLOGY Double-blind, placebo-controlled study with 60 patients with treatment-resistant depression (right-sided cTBS, left-sided iTBS, combined cTBS+iTBS, sham: 15 subjects per group).

PROTOCOL

Target area	Target Area: left/right DLPFC defined as navigated BA9/BA46
Stim. frequency	Stim. frequency: 3-pulse 50 Hz bursts applied at frequency of 5 Hz: (1) right-sided cTBS (1 train of 120 s), (2) left-sided iTBS (trains of 2 s with 10 s pause), (3) combination of both
Stim. intensity	80 % RMT
Number of pulses per session	Number of pulses per session: 1,800 pulses on each site
Treatment	daily, 5 times per week
Number of sessions	10
Assessment	HDRS-17



RESULTS

	Right-sided cTBS	Left-sided iTBS	Combined	Sham
Rate of change of HDRS scores	-22.5 %	-42.3 %	-52.5 %	-17.4 %
Responders	25 %	40 %	66.7 %	13.3 %
Patients with moderate treatment resistance				
Change in HDRS scores	-9.8	-8.6	-15.0	-4.0
Patients with high treatment resistance				
Change in HDRS scores	-1.4	-8.4	-10.6	0.1

CONCLUSIONS

Active theta-burst stimulation is safe, well-tolerated and effective in the treatment of therapy-resistant depression. The number of sessions needed for reduction in HDRS scores seems to be higher in case of greater treatment resistance. Left-sided iTBS and combined cTBS + iTBS are more effective than cTBS in highly refractory depression.

Class of the study: II


Obsessive compulsive disorder (OCD)

Repetitive transcranial magnetic stimulation in the treatment of obsessive-compulsive disorders: Double blind randomized clinical trial

AUTHORS	Elbeh KA, Elserogy YM, Khalifa HE, Ahmed MA, Hafez MH, Khedr EM
SOURCE	Psychiatry Res 238 (2016) 264–269
OBJECTIVE	Effectiveness of 1 Hz and 10 Hz rTMS compared to sham in patients with OCD.
METHODOLOGY	Double-blind, placebo-controlled study with 45 patients (1 Hz, 10 Hz, sham: 15 subjects per group).

PROTOCOL

Target area	right DLPFC (5 cm rule)
Stim. frequency	1 Hz / 10 Hz
Stim. intensity	100 % RMT
Number of pulses per session	2,000
Treatment	daily, 5 times per week
Number of sessions	10 sessions
Assessment	Y-BOCS, HAM-A, CGI-S



RESULTS

Effects on scores of:		1 Hz	10 Hz	sham
Y-BOCS	post-treatment	-12.0*	-6.8*	-1.4***
	3 month follow-up	-11.0*	-2.5***	-2.1***
HAM-A	post-treatment	-12.3*	-7.3*	-1.6***
	3 month follow-up	-11.9*	-3.6***	-3.0***
CGI-S	post-treatment	-2.2*	-1.1**	-0.4***
	3 month follow-up	-2.0*	-0.5***	-0.5***

*p<0.01 **p<0.05 ***p>0.05

CONCLUSIONS

1 Hz rTMS over the right DLPFC shows significant reduction of OCD symptoms after 10 sessions. The improvement is evident in obsessions and compulsions as well as in anxiety. The effect was maintained in the follow-up after 3 months. In conclusion, low frequency rTMS over the right DLPFC is a promising tool for the treatment of OCD.

Class of the study: II


Anxiety

Repetitive transcranial magnetic stimulation of the right dorsal lateral prefrontal cortex in the treatment of generalized anxiety disorder: A randomized, double-blind sham controlled clinical trial

AUTHORS Dilkov D, Hawken ER, Kaludiev E, Milev R
SOURCE Prog Neuropsychopharmacol Biol Psychiatry 78 (2017) 61–65
OBJECTIVE Effectiveness of 20 Hz rTMS compared to sham in patients with generalized anxiety disorder (GAD)
METHODOLOGY Single-blind, placebo-controlled study with 40 patients (15 active rTMS, 25 sham stimulation)

PROTOCOL

Target area	right DLPFC (5 cm rule)
Stim. frequency	20 Hz
Stim. intensity	110 % RMT
Number of pulses per session	3,600
Treatment	5 sessions in week 1 - 4, 3 session in week 5, 2 sessions in week 6
Number of sessions	25
Assessment	HARS, HDRS-21, CGI



RESULTS

Reduction in scores		Active	Sham	Significance
HARS	post-treatment	-25	-1	p<0.001
	6 week follow-up	-26	-1	
HDRS-21	post-treatment	-11	0	p<0.001
	6 week follow-up	-11	1	
CGI	post-treatment	-2	0	p<0.001
	6 week follow-up	-3	0	
Responders*	-1.4	100 %	8 %	

*HARS score reduction >= 50 %

CONCLUSIONS

High-frequency stimulation of the right DLPFC showed clinically significant improvement in anxiety and depression scores in patients with GAD. Follow-up data demonstrate a sustained effect of the response. Considering 50 % of patients do not have a satisfactory response to pharmacotherapies, the treatment successes with rTMS suggest this method is a viable alternative for patients with GAD.


Post-traumatic stress disorder (PTSD)

Unilateral right and bilateral dorsolateral prefrontal cortex transcranial magnetic stimulation in treatment post-traumatic stress disorder: A randomized controlled study

AUTHORS	Ahmadizadeh MJ, Rezaei M
SOURCE	Brain Res Bull 140 (2018) 334–340
OBJECTIVE	Efficacy of bilateral rTMS and unilateral right rTMS compared to sham in patients with PTSD
METHODOLOGY	Single-blind, placebo-controlled study with 65 patients (21 unilateral, 22 bilateral, 22 sham)

PROTOCOL

Target area	bilateral/right DLPFC (Beam F3)
Stim. frequency	20 Hz
Stim. intensity	100 % of RMT
Number of pulses per session	2,400
Treatment	3 times per week for 2 weeks, twice per week for another 2 weeks
Number of sessions	10
Assessment	PCL-M



RESULTS

Change in PCL-M scores	Unilateral	Bilateral	Sham
...after 5 sessions	-2.36	-11.16	-1.8
...after 10 sessions	-21.16	-25.45	-3.62
Responders	41 %	62 %	

*HARS score reduction \geq 50 %

CONCLUSIONS

Bilateral and unilateral right high-frequency rTMS leads to significant reductions of PTSD symptoms. The improvement occurs faster with bilateral stimulation. However, there is no significant difference after 10 sessions. The treatment is safe and well tolerated.

Class of the study: II

Malignant neuropathic pain

Repetitive transcranial magnetic stimulation in neuropathic pain secondary to malignancy: a randomized clinical trial

AUTHORS Khedr EM, Kotb HI, Mostafa MG, Mohamad MF, Amr SA, Ahmed MA, Karim AA, Kamal SMM


SOURCE Eur J Pain 19 (2015) 519–527

OBJECTIVE Efficacy of rTMS over the primary motor cortex (M1) in patients with malignant neuropathic pain.

METHODOLOGY Randomized, placebo-controlled trial including 34 patients (17 active, 17 sham) suffering from malignant neuropathic pain.

PROTOCOL

Target area	hand M1 contralateral to pain
Stim. frequency	20 Hz
Stim. intensity	80 % of RMT
Number of pulses per session	2,000
Treatment	daily, 5 times per week
Number of sessions	10
Assessment	VDS, VAS, LANSS, HAM-D



RESULTS

	Active	Sham
VDS, VAS	Significant difference in pain relief in favor of active rTMS after 10 sessions and 15 days follow-up	
LANSS, HAM-D	Significant difference in pain relief in favor of active rTMS after 10 sessions and 15 days and 1 month follow-up	
Response rates*:		
After sessions	86.6 %	6.6 %
Follow-up 15 days	80 %	6.6 %
Follow-up 1 month	26.6 %	6.6 %

*>30% pain relief

CONCLUSIONS

The results demonstrate that 10 rTMS sessions over the M1 can induce pain relief in malignant neuropathic pain. A high proportion of patients in the active group responded to HF rTMS with lasting effect up to 15 days. In about one quarter of patients the effect remained stable 1 month after treatment.

Class of the study: II


Chronic neuropathic pain

Time course of the response to navigated repetitive transcranial magnetic stimulation at 10 Hz in chronic neuropathic pain

AUTHORS	Lawson McLean A, Frank S, Zafar N, Waschke A, Kalff R, Reichart R
SOURCE	Neurol Res 40 (2018) 564–572
OBJECTIVE	Evaluation of the time to response and effectiveness of navigated rTMS at a frequency of 10 Hz in patients with chronic neuropathic pain.
METHODOLOGY	Prospective study including 48 patients with unilateral chronic neuropathic pain.

PROTOCOL

Target area	contralateral motor cortex (atypical facial, upper limb and torso pain: hand area; lower limb pain: lower limb area)
Stim. frequency	10 Hz
Stim. intensity	80 % of RMT
Number of pulses per session	1,000
Treatment	daily, 4 - 5 times per week
Number of sessions	9
Assessment	VAS, German Pain Questionnaire



RESULTS

Response rate	58.3 %
Responders in 6 weeks follow-up	41.6 %
VAS	
Post stimulation 6 weeks follow-up	significant pain reduction in the responder group (p<0.001)
Pain reduction in the responder group:	
Atypical facial pain	54.8 %
Central post-stroke pain	66.3 %
Neuropathic lower limb pain	42.9 %
others	35.5 %
Pain relief onset	after 3 - 41 days

CONCLUSIONS

Clinical benefit of navigated rTMS in treatment of chronic pain has been demonstrated. In particular, patients with chronic painful post-traumatic trigeminal neuropathy responded well to high-frequency stimulation. Patients with a mean pain history of less than five years benefited significantly from this treatment, so early treatment with repetitive TMS should be encouraged.


Motor disorders - stroke

Daily repetitive transcranial magnetic stimulation for poststroke upper limb paresis in the subacute period

AUTHORS	Hosomi K, Morris S, Sakamoto T, Taguchi J, Maruo T, Kageyama Y, Kinoshita Y, Goto Y, Shimokawa T, Komaya T, Saitoh Y
SOURCE	J Stroke Cerebrovasc Dis 25 (2016) 1655–1664
OBJECTIVE	Efficacy and safety of daily rTMS in motor recovery of subacute stroke patients.
METHODOLOGY	Double-blind, placebo-controlled study with 41 patients (20 active, 21 sham). rTMS accompanied a regular rehabilitation program. Daily rehabilitation consisted of 8 sessions a 20 min including physical and occupational therapy.

PROTOCOL

Target area	Hand M1 (C3)
Stim. frequency	5 Hz
Stim. intensity	90 % RMT
Number of pulses per session	500
Treatment	daily, 5 times per week
Number of sessions	10
Assessment	BS, FMA, Handgrip strength, NIHSS, FIM



RESULTS

	Active	Sham
BS	Significant improvement in hand score, Significant difference between active and sham for hand score, Significant improvement in arm score, Improvement in lower limb score	Improvement in arm and lower limb score
FMA, NIHSS, FM	Significant improvement compared to baseline	
Hand grip strength	Significant improvement No significant group difference	Improvement
Adverse events	No serious adverse events	

CONCLUSIONS

Daily HF rTMS of the ipsilesional M1 is tolerable and facilitates motor recovery in the paralytic hand of subacute stroke patients. The effectiveness of physical and occupational therapy in the framework of an intense rehabilitation program in the subacute stage after stroke can be increased by accompanying rTMS.

Class of the study: II

Motor disorders - stroke

Effect of combined low-frequency repetitive transcranial magnetic stimulation and virtual reality training on upper limb function in subacute stroke: a double-blind randomized controlled trial

AUTHORS Zheng C, Liao W, Xia W


SOURCE J Huazhong Univ Sci Technolog Med Sci 35 (2015) 248–254

OBJECTIVE Effect of combined low-frequency rTMS and virtual reality training in stroke patients.

METHODOLOGY Double-blind, placebo-controlled trial including 108 patients with hemiplegia after stroke (55 active, 53 sham) on average 19 days poststroke. The rTMS treatment was combined with a rehabilitation program which included VR training, occupational therapy and physiotherapy.

PROTOCOL

Target area	Hand M1 (C3)
Stim. frequency	1 Hz
Stim. intensity	90% RMT
Number of pulses per session	1,800
Treatment	daily, 6 times per week
Number of sessions	24
Assessment	U-FMA, WMFT, MBI, SF-36



RESULTS

	Experimental group			Control group		
	Before treatment	Four weeks after treatment	dif	Before treatment	Four weeks after treatment	dif
WMFT scores	32.4	51.8	25.7	31.6	44.7	13.1
MBI scores	52.6	87.2	34.6	53.4	71.6	18.2
SF-36 scores (PF)	34.4	65.6	31.2	33.6	47.4	13.8

CONCLUSIONS

LF rTMS combined with intensive rehabilitation programs (physiotherapy, occupational therapy, VR training) can effectively improve the upper limb function, the living activity, and the quality of life in patients with hemiplegia following stroke. The VR training seems to be particularly boosted by LF rTMS finished only a few minutes before starting the VR training.

Class of the study: I


Aphasia - stroke

The persistent and broadly modulating effect of inhibitory rTMS in nonfluent aphasic patients: a sham-controlled, double-blind study

AUTHORS	Tsai PY, Wang CP, Ko JS, Chung YM, Chang YW, Wang JX
SOURCE	Neurorehabil Neural Repair 28 (2014) 779–787
OBJECTIVE	Effectiveness of rTMS on language performance in stroke patients with nonfluent aphasia. Identification of characteristics of patients predisposed to benefit most.
METHODOLOGY	Randomized, placebo-controlled study with 56 patients with nonfluent aphasia 18 months poststroke on average (33 LF rTMS, 23 sham) combined with 1 hour of speech therapy.

PROTOCOL

Target area	contralateral pars triangularis
Stim. frequency	1 Hz
Stim. intensity	90 % RMT
Number of pulses per session	600
Treatment	daily, 5 times per week
Number of sessions	10
Assessment	CCAT, Picture naming test



RESULTS

Overall CCAT score after 10 sessions	Significant improvement compared to baseline, $p < 0.001$ Significant improvement compared to sham, $p < 0.001$
Subscores of CCAT after 10 sessions	Significant improvement compared to sham in conversation, $p < 0.032$, description, $p < 0.024$, expression, $p < 0.002$, repetition, $p < 0.023$
Overall CCAT score in follow-up 3 month	Significant improvement compared to baseline, $p < 0.008$
Subscores of CCAT in follow-up 3 month	Significant improvement compared to baseline in description, $p < 0.031$ expression, $p < 0.026$, repetition, $p < 0.013$
Impact on therapeutic outcome	<ul style="list-style-type: none"> - outcome independent of severity of non-fluent aphasia - outcome independent of time post-stroke - history of diabetes mellitus: negative impact on improvement, $p < 0.024$ - high contralateral corticomotor excitability: better outcome of inhibitory rTMS, $p < 0.006$

CONCLUSIONS

Downregulating the circuitry of the right pars triangularis (Ptr) through inhibitory rTMS, achieves a persistent and broadly modulating effect, irrespective of aphasia severity and subtype. Patients with lower RMT in the right motor system seem to benefit the most.

Class of the study: II

Multiple sclerosis

High-frequency repetitive transcranial magnetic stimulation and intermittent thetaburst stimulation for spasticity management in secondary progressive multiple sclerosis

AUTHORS Korzhova J, Bakulin I, Sinitsyn D, Poydasheva A, Suponeva N, Zakharova M, Piradov M

SOURCE Eur J Neurol 26 (2019) 680-686


OBJECTIVE Comparison of HF-rTMS (20 Hz) and iTBS on the level of spasticity and concomitant symptoms in patients with secondary progressive multiple sclerosis.

METHODOLOGY Randomized, placebo-controlled study with 34 patients (12 HF-rTMS, 12 iTBS, 10 sham, each combined with physical therapy).

PROTOCOL

Target area	Bilateral leg M1
Stim. frequency	20 Hz / iTBS*
Stim. intensity	80 % MSO**
Number of pulses per session	1,600 / 1,200
Treatment	daily, 5 times per week
Number of sessions	10 sessions
Assessment	MAS, NAS, SESS, MFIS

* bursts at 5Hz containing 3 pulses at 35 Hz **Maximum stimulant power



RESULTS

	20 Hz	iTBS	Sham
MAS, SESS, NAS after 10 sessions	Significant reduction of spasticity in MAS, NAS, not significant in SESS	Significant reduction of spasticity	No significant reduction of spasticity
NAS, SESS follow-up 2 weeks	Significant reduction of spasticity	Significant reduction of spasticity	No significant reduction of spasticity
NAS, SESS follow-up 12 weeks	No significant reduction of spasticity	Significant reduction of spasticity	No significant reduction of spasticity
Pain level (spasticity related)	Significant reduction of pain at end of sessions and in follow-up 2 weeks, not significant in follow-up 12 weeks	No significant reduction	No significant reduction
Fatigue (MFIS)	Significant reduction of fatigue at end of sessions and in follow-up 2 weeks, gradual increase until follow-up 12 weeks	No significant reduction	No significant reduction

CONCLUSIONS

The results show that HF-rTMS and iTBS significantly reduce spasticity, in contrast to sham stimulation. Some evidence was found in favor of a longer-lasting effect of iTBS. A reduction in pain and fatigue could be achieved by HF-rTMS but did not result from iTBS.

Class of the study: II



neuroCademy 

Would you like to learn more about TMS?

Training Courses by neuroCademy help therapists deliver TMS in their own clinical setting. The courses follow the neuroCare Group's proven and highly effective protocols and feature specialized knowledge on localizing the DLPFC and enhancing treatment with Psychotherapy. Training is suited to Psychiatrists, Psychologists and Mental Health Nurses.

www.neuroCademy.com

Training & Science

NEW: E-Learning at neuroCademy

neuroCademy courses about TMS and other neuromodulation methods are now offered online as well. Attend our self-placed online courses on the neuroCademy Learning Management System "LMS". The courses are modular and include on-demand videos as well as moderated sessions.

<https://lms.neuroCademy.com>

Information, advice and registration:

Europe: +49 3677 68 979 0

Australia: +61 2 8317 5032

neurocademy@neurocaregroup.com

Abbreviations

BS	Brunnstrom stages
CCAT	Concise Chinese Aphasia Test
CGI-S	Clinical Global Impression-Severity scale
DASS	Depressions-Angst-Stress-Skala
DLPFC	dorsolateral prefrontal cortex
FIM	Functional Independence Measure
FMA	Fugl-Meyer Assessment
GAD	general anxiety disorder
HAM-A	Hamilton Anxiety Rating Scale
HAM-D	Hamilton Rating Scale for Depression
HF-rTMS	high-frequency rTMS
iTBS	intermittent theta-burst stimulation
LANSS	Leeds Assessment of Neuropathic Symptoms and Signs
MADR	Montgomery–Asberg Depression Rating Scale
MAS	Modified Ashworth Scale
MBI	Modified Barthel Index
MFIS	Modified Fatigue Impact Scale
MSO	Maximum Stimulant Output
NAS	Numerical Analog Scale
NIHSS	National Institutes of Health Stroke Scale
OCD	obsessive compulsive disorder
SESS	Subjective Evaluating Spasticity Scale
SF-36	Short Form Health Survey Questionnaire
TPC	temporo-polar cortex
U-FMA	Fugl-Meyer Assessment for upper limb
VAS	Visual Analogue Scale
VDS	Verbal Descriptor Scale
WMFT	Wolf Motor Function Test
Y-BOCS	Yale-Brown obsessive compulsive scale

Imprint

neurocare group AG · Albert-Einstein-Straße 3 · 98693 Ilmenau · Germany
Administrative office: Rindermarkt 7 · 80331 Munich · Germany
Tel: +49 (0)89 3564 767 0 · Fax: +49 (0)89 2154 712.99 1
eMail: info@neurocaregroup.com · Web: www.neurocaregroup.com

Images: neurocare group AG

Copyright: neurocare group AG 2021

The use or publication of contained texts or pictures is strictly prohibited.
Exceptions require the written approval of neurocare group AG.

neuroConn 
Technology

neuroCademy 
Training & Science

neuroCare 
Clinics

neurocare group AG

Albert-Einstein-Straße 3 | 98693 Ilmenau | Germany

phone: +49 3677 68 979 0

e-mail: sales@neurocaregroup.com

web: www.neurocaregroup.com

administrative office: Rindermarkt 7 | 80331 Munich | Germany

phone: +49 89 3564 767 0