


# Transcranial Magnetic Stimulation TMS

Effectiveness in studies





# Content

- 4** rTMS - Overview and evidence-based guidelines
- 6** Depression
- 12** Obsessive Compulsive Disorder (OCD)
- 14** Anxiety
- 16** Post-traumatic stress disorder (PTSD)
- 18** Pain
- 22** Stroke - motor disorders
- 26** Stroke - Aphasia
- 28** Multiple Sclerosis
  
- 30** Abbreviations
- 31** Training & Science

## **Imprint**

neurocare group AG

Albert-Einstein-Str. 3, 98693 Ilmenau, Germany

Phone: +48 (3677) 68 979 0

E-mail: [info@neurocaregroup.com](mailto:info@neurocaregroup.com) • web: [www.neurocaregroup.com](http://www.neurocaregroup.com)

Administrative Office: Rindermarkt 7 • 80331 Munich • Germany

Images: neurocare group AG, Copyright: neurocare group AG 2020

The use or publication of contained text or images is strictly prohibited.

Exceptions require the written approval of neurocare group AG.



# TMS: Overview – evidence-based guidelines

Evidence-based guidelines on the therapeutic use of repetitive transcranial magnetic stimulation (rTMS): An update (2014–2018).

## Authors

Lefaucheur JP, Aleman A, Baeken C, Benninger DH, Brunelin J, Di Lazzaro V, Filipovic SR, Grefkes C, Hasan A, Hummel FC, Jääskeläinen SK, Langguth B, Leocani L, Londero A, Nardone R, Nguyen JP, Nyffeler T, Oliveira-Maia AJ, Oliviero A, Padberg F, Palm U, Paulus W, Poulet E, Quartarone A, Rachid F, Rektorová I, Rossi S, Sahlsten H, Schecklmann M, Szekely D, Ziemann U

## Source

Clinical Neurophysiology, 131 (2020), pp 474–528

## Objective

overview of the latest research results on the therapeutic use of rTMS and evaluation of the scientific evidence.

## Methodology

classification of studies in terms of criteria randomized/controlled, number of patients, clear statements on procedure, primary outcome, exclusion criteria, statistical analysis, patient groups. Summary of outcomes

## Results

---

<b>Level A</b>	depression (high-frequency rTMS) • neuropathic pain • post-stroke contralateral hand motor recovery in post-acute stage
<b>Level B</b>	depression (low-frequency rTMS, bilateral rTMS • post-traumatic stress disorder • fibromyalgia • Parkinson's disease: depressive symptoms • Parkinson's disease: motor symptoms • post-stroke ipsilateral hand motor recovery in post-acute stage • post-stroke nonfluent aphasia in chronic stage • lower limb spasticity in multiple sclerosis
<b>Level C</b>	obsessive compulsive disorder • addiction/craving • schizophrenia: auditory hallucinations & negative symptoms • tinnitus • hemispatial neglect in post-acute stage after stroke • pilepsy • Alzheimer's disease: cognitive function • complex regional pain syndrome – type I

---

## Conclusions

A review of hundreds of scientific studies has shown that rTMS can provide significant improvement in various neurological and psychiatric disorders. rTMS is scientifically proven to be definitely effective in the treatment of neuropathic pain and depression as well as in the hand motor recovery from stroke. Furthermore, scientific studies show that rTMS is probably effective for patients with motor impairment, Parkinson's disease, multiple sclerosis, post-traumatic stress disorder or aphasia.

# Depression

Simultaneous rTMS and psychotherapy in major depressive disorder: Clinical outcomes and predictors from a large naturalistic study.

## Authors

Donse L, Padberg F, Sack AT, Rush AJ, Arns M

## Source


Brain Stimulation 11 (2018) 337-345

## Objective

evaluation of feasibility and clinical outcome of rTMS in combination with psychotherapy in major depressive disorder. Identification of clinical predictors of response and remission

## Methodology

naturalistic study including 196 patients with depression, rTMS treatment with simultaneous psychotherapy. 97% had been unsuccessfully treated with medication at least once before.



## Protocol

<b>target area</b>	left DLPFC / right DLPFC
<b>stim. frequency</b>	10 Hz / 1 Hz
<b>stim. intensity</b>	110 - 120% RMT
<b>number of pulses per session</b>	1,500 / 1,200
<b>treatment</b>	2 - 10 times per week
<b>scope</b>	21 treatments on average
<b>assessment</b>	BDI, DASS, Reaktion defined as $\geq 50\%$ reduction in BDI score



## Results

---

### assessment

---

**response rate** 66,3%

---

**remission** 56,0%

---

**BDI** baseline 31.3 / Last session 14.1 / 55.9% reduction,  $p < 0.001$

---

**DASS D** baseline 28.6 / last session 12.1

**DASS A** baseline 13.7 / last session 5.6

**DASS S** baseline 22.4 / last session 10.1

---

**follow-up 6 months** group of responders

---

**BDI** baseline 29.4 / last session 8.0 / follow-up 13.8

**DASS depression** baseline 27.0 / last session 6.5 / follow-up 11.8

**DASS anxiety** baseline 14.1 / last session 4.0 / follow-up 7.2

**DASS stress** baseline 22.3 / last session 7.3 / follow-up 12.5

---

## Conclusions

The combined therapy rTMS + psychotherapy led to high response and remission rates. In addition, the effect also showed good stability in the follow-up. Considering the high proportion of treatment resistance against antidepressants in the sample, the result is of clinical relevance.

# Depression

Efficacy and Safety of Transcranial Magnetic Stimulation in the Acute Treatment of Major Depression: A Multisite Randomized Controlled Trial.

## Authors

O'Reardon JP, Solvason HB, Janicak PG, Sampson S, Isenberg KE, Nahas Z, McDonald WM, Avery D, Fitzgerald PB, Loo C, Demitrack MA, George MS, Sackheim HA

## Source


Brain Stimulation 11 (2018) 337-345

## Objective

effectiveness and safety of TMS over left DLPFC in the acute treatment of major depression

## Methodology

double-blind, placebo-controlled study with 301 patients with depression (155 active rTMS, 146 placebo stimulation). Patients were medication-free and had shown medication resistance prior to rTMS treatment.



## Protocol

<b>target area</b>	left DLPFC using 5 cm rule
<b>stim. frequency</b>	10 Hz
<b>stim. intensity</b>	120 % RMT
<b>number of pulses per session</b>	3,000
<b>treatment</b>	every day, 5 times per week
<b>scope</b>	max. 30 treatments
<b>assessment</b>	MADRS, HAMD17, HAMD24



## Results

assessment	significance active vs. placebo	active stimulation	placebo stimulation
MADRS response rate*	<b>P &lt; 0.01</b>	23.9 %	12.3 %
HAMD24 response rate	<b>P &lt; 0.05</b>	23.9 %	15.1 %
MADRS remission rate**	<b>P &lt; 0.05</b>	14.2 %	5.5 %
HAMD24 remission rate***	<b>P &lt; 0.05</b>	17.4 %	8.2 %
<b>adverse events</b>			
scalp discomfort/ pain reported as mild to moderate		35.8 %; diminished after first 5 sessions	3.8 %
increase in suicidality		<b>1 event</b>	<b>10 events</b>

\* 50% improvement from baseline.

\*\* total score <10

\*\*\* total score <8

\*\*\*\* total score <11

## Conclusions

TMS administered over the left DLPFC using the parameters reported here for a period of up to 6 weeks was effective in treating major depression with a good tolerability profile. These results indicate that rTMS offers clinicians a novel alternative in the treatment of this disorder.

# Depression

Efficacy of prefrontal theta-burst stimulation in refractory depression: a randomized placebo-controlled study.

## Authors

Li CT, Chen MH, Juan CH, Huang HH, Chen LF, Hsieh JC, Tu PC, Bai YM, Tsai SJ, Lee YC, Su TP

## Source


Brain 137 (2014) 2088–2098

## Objective

antidepressant effect of right pre-frontal continuous, left-prefrontal intermittent and combined theta-burst stimulation compared to placebo group of patients showing different levels of antidepressant treatment resistance

## Methodology

double-blind, placebo-controlled study with 60 patients with treatment-resistant depression (right-sided cTBS, left-sided iTBS, combined cTBS+iTBS, placebo: 15 subjects per group)



## Protocol

**target area** left/right DLPFC defined as navigated BA9/BA46

**stim. frequency** 3-pulse 50 Hz bursts applied at frequency of 5 Hz:  
 (1) right-sided cTBS (1 train of 120 s),  
 (2) left-sided iTBS (trains of 2 s with 10 s pause),  
 (3) combination of both

**stim. intensity** 80% RMT

**number of pulses per session** 1,800 pulses per side

**treatment** every day, 5 times per week

**scope** 10 treatments

**assessment** HDRS-17

## Results

	right-sided cTBS	left-sided iTBS	combination	placebo stimulation
rate of change of HDRS scores	-22.5%	-42.3%	-52.5%	-17.4%
responders	25%	40%	66.7%	13.3%
<b>patients with moderate treatment resistance</b>				
change in HDRS scores	-9.8	-8.6	-15.0	-4.0
<b>patients with high treatment resistance</b>				
change in HDRS scores	-1.4	-8.4	-10.6	0.1

## Conclusions

Active theta-burst stimulation is safe, well-tolerated and effective in the treatment of therapy-resistant depression. The number of sessions needed for reduction in HDRS scores seems to be higher in case of greater treatment resistance. Left-sided iTBS and combined cTBS + iTBS are more effective than cTBS in highly refractory depression.

# Obsessive-compulsive disorder (OCD)

Repetitive transcranial magnetic stimulation in the treatment of obsessive-compulsive disorders: Double blind randomized clinical trial.

## Authors

Elbeh KA, Elserogy YM, Khalifa HE, Ahmed MA, Hafez MH, Khedr EM

## Source


Psychiatry Res 238 (2016) 264–269

## Objective

effectiveness of 1 Hz and 10 Hz rTMS compared to placebo in patients with OCD

## Methodology

double-blind, placebo-controlled study with 45 patients (1 Hz, 10 Hz, placebo: 15 subjects per group)



## Protocol

<b>target area</b>	right DLPFC (5 cm rule)
<b>stim. frequency</b>	1 Hz / 10 Hz
<b>stim. intensity</b>	100 % RMT
<b>number of pulses per session</b>	2,000
<b>treatment</b>	every day, 5 time per week
<b>scope</b>	10 treatments
<b>assessment</b>	Y-BOCS, HAM-A, CGI-S

## Results

		1 Hz	10 Hz	placebo stimulation
Y-BOCS	post-treatment	-12.0*	-6.8*	-1.4***
	3 month follow-up	-11.0*	-2.5***	-2.1***
HAM-A	post-treatment	-12.3*	-7.3*	-1.6***
	3 month follow-up	-11.9*	-3.6***	-3.0***
CGI-S	post-treatment	-2.2*	-1.1**	-0.4***
	3 month follow-up	-2.0*	-0.5***	-0.5***

\* p<0.01

\*\* p<0.05

\*\*\* p>0.05

## Conclusions

1 Hz rTMS over the right DLPFC shows significant reduction of OCD symptoms after 10 sessions. The improvement is evident in obsessions and compulsions as well as in anxiety. The effect was maintained in the follow-up after 3 months. In conclusion, low frequency rTMS over the right DLPFC is a promising tool for the treatment of OCD.

# Anxiety

Repetitive transcranial magnetic stimulation of the right dorsal lateral prefrontal cortex in the treatment of generalized anxiety disorder: A randomized, double-blind placebo controlled clinical trial.

## Authors

Dilkov D, Hawken ER, Kaludiev E, Milev R

## Source


Prog Neuropsychopharmacol Biol Psychiatry 78 (2017) 61–65

## Objective

effectiveness of 20 Hz rTMS compared to placebo in patients with generalized anxiety disorder (GAD)

## Methodology

single-blind, placebo-controlled study with 40 patients (15 active rTMS, 25 placebo stimulation)



## Protocol

<b>target area</b>	right DLPFC (5 cm rule)
<b>stim. frequency</b>	20 Hz
<b>stim. intensity</b>	110 % RMT
<b>number of pulses per session</b>	3,600
<b>treatment</b>	5 treatments in week 1 - 4 3 treatments in week 5 2 treatments in week 6
<b>scope</b>	25 treatments
<b>assessment</b>	HARS, HDRS-21, CGI

## Results

reductions in scores		active stimulation	placebo stimulation	significance
<b>HARS</b>	post-treatment	-25	-1	<b>p&lt;0.001</b>
	6 week follow-up	-26	-1	
<b>HDRS-21</b>	post-treatment	-11	0	<b>p&lt;0.001</b>
	6 week follow-up	-11	1	
<b>CGI</b>	post-treatment	-2	0	<b>p&lt;0.001</b>
	6 week follow-up	-3	0	
<b>responders*</b>	-1.4	<b>100%</b>	8%	

\* HARS score reduction >= 50 %

## Conclusions

High-frequency stimulation of the right DLPFC showed clinically significant improvement in anxiety and depression scores in patients with GAD. Follow-up data demonstrate a sustained effect of the response. Considering 50 % of patients do not have a satisfactory response to pharmacotherapies, the treatment successes with rTMS suggest this method is a viable alternative for patients with GAD.



# Post-traumatic stress disorder (PTBS)

Unilateral right and bilateral dorsolateral prefrontal cortex transcranial magnetic stimulation in treatment post-traumatic stress disorder: A randomized controlled study.

## Authors

Ahmadizadeh MJ, Rezaei M

## Source


Brain Res Bull 140 (2018) 334–340

## Objective

efficacy of bilateral rTMS and unilateral right rTMS compared to placebo in patients with PTSD

## Methodology

single-blind, placebo-controlled study with 65 patients (21 unilateral, 22 bilateral, 22 placebo)



## Protocol

<b>target area</b>	bilateral / right DLPFC (Beam F3)
<b>stim. frequency</b>	20 Hz
<b>stim. intensity</b>	100 % of RMT
<b>number of pulses per session</b>	2,400
<b>treatment</b>	3 treatments per week for 2 weeks, twice per week for another 2 weeks
<b>scope</b>	10 treatments
<b>assessment</b>	PCL-M

## Results

change in PCL-M scores	unilateral	bilateral	placebo stimulation
...after 5 sessions	-2.36	-11.16	-1.8
...after 10 sessions	-21.16	-25.45	-3.62
responders*	41 %	62 %	

## Conclusions

Bilateral and unilateral right high-frequency rTMS leads to significant reductions of PTSD symptoms. The improvement occurs faster with bilateral stimulation. However, there is no significant difference after 10 sessions. The treatment is safe and well tolerated.

# Malignant neuropathic pain

Repetitive transcranial magnetic stimulation in neuropathic pain secondary to malignancy: a randomized clinical trial.

## Authors

Khedr EM, Kotb HI, Mostafa MG, Mohamad MF, Amr SA, Ahmed MA, Karim AA, Kamal SMM

## Source


Eur J Pain 19 (2015) 519–527

## Objective

efficacy of rTMS over the primary motor cortex (M1) in patients with malignant neuropathic pain

## Methodology

randomized, placebo-controlled trial including 34 patients (17 active, 17 placebo) suffering from malignant neuropathic pain



## Protocol

<b>target area</b>	hand: M1 contralateral to side of pain
<b>stim. frequency</b>	20 Hz
<b>stim. intensity</b>	80% RMT
<b>number of pulses per session</b>	2,000
<b>treatment</b>	every day, 5 times per week
<b>scope</b>	10 treatments
<b>assessment</b>	VDS, VAS, LANSS, HAM-D

## Results

	active stimulation	placebo stimulation
<b>VDS, VAS</b>	significant difference in pain relief in favor of active rTMS after 10 sessions and 15 days follow-up	
<b>LANSS, HAM-D</b>	significant difference in pain relief in favor of active rTMS after 10 sessions and 15 days and 1 month follow-up	
<b>response rate*:</b>		
<b>after sessions</b>	<b>86.6 %</b>	6.6 %
<b>follow-up 15 days</b>	<b>80 %</b>	6.6 %
<b>follow-up 1 month</b>	<b>26.6 %</b>	6.6 %

\* > 30% pain relief

## Conclusions

The results demonstrate that 10 rTMS sessions over the M1 can induce pain relief in malignant neuropathic pain. A high proportion of patients in the active group responded to HF rTMS with lasting effect up to 15 days. In about one quarter of patients the effect remained stable 1 month after treatment.

# Chronic neuropathic pain

Time course of the response to navigated repetitive transcranial magnetic stimulation at 10 Hz in chronic neuropathic pain

## Authors

Lawson McLean A, Frank S, Zafar N, Waschke A, Kalff R, Reichart R

## Source


Neurol Res 40 (2018) 564–572

## Objective

evaluation of the time to response and effectiveness of navigated rTMS at a frequency of 10 Hz in patients with chronic neuropathic pain

## Methodology

prospective study including 48 patients with unilateral chronic neuropathic pain



## Protocol

<b>target area</b>	contralateral motor cortex (atypical facial, upper limb and torso pain: hand area; lower limb pain: lower limb area)
<b>stim. frequency</b>	10 Hz
<b>stim. intensity</b>	80% of RMT
<b>number of pulses per session</b>	1,000
<b>treatment</b>	daily, 4 - 5 times per week
<b>scope</b>	9 treatments
<b>assessment</b>	VAS, German Pain Questionnaire

## Results

---

<b>response rate</b>	<b>58.3%</b>
<b>responders in 6 weeks follow-up</b>	<b>41.6%</b>

---

<b>VAS</b>	<b>significant pain reduction in the responder group</b>
<b>post stimulation</b>	<b>(p&lt;0.001)</b>
<b>6 weeks follow-up</b>	

---

<b>pain reduction in the responder group:</b>	
<b>atypical facial pain</b>	54.8%
<b>central post-stroke pain</b>	66.3%
<b>neuropathic lower limb pain</b>	42.9%
<b>others</b>	35.5%

---

<b>pain relief onset</b>	after 3 - 41 days
--------------------------	-------------------

---

## Conclusions

Clinical benefit of navigated rTMS in treatment of chronic pain has been demonstrated. In particular, patients with chronic painful post-traumatic trigeminal neuropathy responded well to high-frequency stimulation. Patients with a mean pain history of less than five years benefited significantly from this treatment, so early treatment with repetitive TMS should be encouraged.

# Motor disorders – stroke

Daily repetitive transcranial magnetic stimulation for poststroke upper limb paresis in the subacute period.

## Authors

Hosomi K, Morris S, Sakamoto T, Taguchi J, Maruo T, Kageyama Y, Kinoshita Y, Goto Y, Shimokawa T, Komaya T, Saitoh Y

## Source


J Stroke Cerebrovasc Dis 25 (2016) 1655–1664

## Objective

efficacy and safety of daily rTMS in motor recovery of subacute stroke patients

## Methodology

double-blind, placebo-controlled study with 41 patients (20 active, 21 placebo); rTMS accompanied a regular rehabilitation program; daily rehabilitation consisted of 8 sessions a 20 min including physical and occupational therapy



## Protocol

<b>target area</b>	hand M1 (C3)
<b>stim. frequency</b>	5 Hz
<b>stim. intensity</b>	90% RMT
<b>number of pulses per session</b>	500
<b>treatment</b>	every day, 5 time per week
<b>scope</b>	10 treatments
<b>assessment</b>	BS, FMA, handgrip strength, NIHSS, FIM



## Results

	<b>active stimulation</b>	<b>placebo stimulation</b>
<b>BS</b>	significant improvement in hand score, significant difference between active and placebo for hand score, significant improvement in arm score, improvement in lower limb score	improvement in arm and lower limb score
<b>FMA, NIHSS, FM</b>	significant improvement compared to baseline	
<b>hand grip strength</b>	significant improvement no significant group difference	improvement
<b>adverse events</b>	no serious adverse events	

## Conclusions

Daily HF rTMS of the ipsilesional M1 is tolerable and facilitates motor recovery in the paralytic hand of subacute stroke patients. The effectiveness of physical and occupational therapy in the framework of an intense rehabilitation program in the subacute stage after stroke can be increased by accompanying rTMS.

# Motor disorders – stroke

Effect of combined low-frequency repetitive transcranial magnetic stimulation and virtual reality training on upper limb function in subacute stroke: a double-blind randomized controlled trial.

## Authors

Zheng C, Liao W, Xia W

## Source


J Huazhong Univ Sci Technolog  
Med Sci 35 (2015) 248–254

## Objective

effect of combined low-frequency rTMS and virtual reality training in stroke patients

## Methodology

double-blind, placebo-controlled trial including 108 patients with hemiplegia after stroke (55 active, 53 placebo) on average 19 days poststroke. The rTMS treatment was combined with a rehabilitation program which included VR training, occupational therapy and physiotherapy



## Protocol

<b>target area</b>	hand M1 (C3)
<b>stim. frequency</b>	1 Hz
<b>stim. intensity</b>	90% RMT
<b>number of pulses per session</b>	1,800
<b>treatment</b>	every day, 6 times per week
<b>scope</b>	24 treatments
<b>assessment</b>	U-FMA, WMFT, MBI, SF-36

## Results

	active stimulation			placebo stimulation		
	before treatment	four weeks after treatment	dif.	before treatment	four weeks after treatment	dif.
<b>WMFT scores</b>	32.4	51.8	<b>25.7</b>	31.6	44.7	<b>13.1</b>
<b>MBI scores</b>	52.6	87.2	<b>34.6</b>	53.4	71.6	<b>18.2</b>
<b>SF-36 scores (PF)</b>	34.4	65.6	<b>31.2</b>	33.6	47.4	<b>13.8</b>

## Conclusions

LF rTMS combined with intensive rehabilitation programs (physiotherapy, occupational therapy, VR training) can effectively improve the upper limb function, the living activity, and the quality of life in patients with hemiplegia following stroke. The VR training seems to be particularly boosted by LF rTMS finished only a few minutes before starting the VR training.

# Aphasia – stroke

The persistent and broadly modulating effect of inhibitory rTMS in nonfluent aphasic patients: a placebo-controlled, double-blind study.

## Authors

Tsai PY, Wang CP, Ko JS, Chung YM, Chang YW, Wang JX

## Source


Neurorehabil Neural Repair 28 (2014) 779–787

## Objective

effectiveness of rTMS on language performance in stroke patients with nonfluent aphasia; identification of characteristics of patients predisposed to benefit most

## Methodology

randomized, placebo-controlled study with 56 patients with nonfluent aphasia 18 months poststroke on average (33 LF rTMS, 23 placebo) combined with 1 hour of speech therapy



## Protocol

<b>target area</b>	contralesional pars triangularis
<b>stim. frequency</b>	1 Hz
<b>stim. intensity</b>	90% RMT
<b>number of pulses per session</b>	600
<b>treatment</b>	every day, 5 times per week
<b>scope</b>	10 treatments
<b>assessment</b>	CCAT, picture naming test

## Results

<b>overall CCAT score after 10 sessions</b>	significant improvement compared to baseline, $p < 0.001$ <b>significant improvement compared to placebo, <math>p &lt; 0.001</math></b>
<b>subscores of CCAT after 10 sessions</b>	<b>significant improvement compared to placebo in conversation</b> in conversation; $p < 0.032$ , description, $p < 0.024$ , expression, $p < 0.002$ , repetition, $p < 0.023$
<b>overall CCAT score in follow-up 3 month</b>	<b>significant improvement compared to baseline, <math>p &lt; 0.008</math></b>
<b>subscores of CCAT in follow-up 3 month</b>	significant improvement compared to baseline in description, $p < 0.031$ expression, $p < 0.026$ , repetition, $p < 0.013$
<b>impact on therapeutic outcome</b>	<b>- outcome independent of severity of non-fluent aphasia</b> <b>- outcome independent of time post-stroke</b> - history of diabetes mellitus: negative impact on improvement, $p < 0.024$ - high contralesional corticomotor excitability: better outcome of inhibitory rTMS, $p < 0.006$

## Conclusions

Down-regulating the circuitry of the right pars triangularis (PT<sub>r</sub>) through inhibitory rTMS, achieves a persistent and broadly modulating effect, irrespective of aphasia severity and subtype. Patients with lower RMT in the right motor system seem to benefit the most.

# Multiple Sclerosis

High-frequency repetitive transcranial magnetic stimulation and intermittent thetaborst stimulation for spasticity management in secondary progressive multiple sclerosis.

## Authors

Korzhoва J, Bakulin I, Sinitsyn D, Poydasheva A, Suponeva N, Zakharova M, Piradov M

## Source


Eur J Neurol 26 (2019) 680-686

## Objective

comparison of HF-rTMS (20 Hz) and iTBS on the level of spasticity and concomitant symptoms in patients with secondary progressive multiple sclerosis

## Methodology

randomized, placebo controlled study with 34 patients (12 HF-rTMS, 12 iTBS, 10 placebo, each combined with physical therapy)



## Protocol

<b>target area</b>	bilateral leg M1
<b>stim. frequency</b>	20 Hz / iTBS*
<b>stim. intensity</b>	80% MSO**
<b>number of pulses per session</b>	1,600 / 1,200
<b>treatment</b>	every day, 5 times per week
<b>scope</b>	10 treatments
<b>assessment</b>	MAS, NAS, SESS, MFIS

\* bursts at 5 Hz containing 3 pulses at 35 Hz

\*\* \*\*\*maximum stimulant power

## Results

	20 Hz	iTBS	placebo stimulation
<b>MAS, SESS, NAS after 10 session</b>	significant reduction of spasticity in MAS, NAS, not significant in SESS	significant reduction of spasticity	no significant reduction of spasticity
<b>NAS, SESS follow-up 2 weeks</b>	significant reduction of spasticity	significant reduction of spasticity	no significant reduction of spasticity
<b>NAS, SESS follow-up 12 week</b>	no significant reduction of spasticity	<b>significant reduction of spasticity</b>	no significant reduction of spasticity
<b>pain level (spasticity related)</b>	<b>significant reduction of pain</b> at end of sessions and in follow-up 2 weeks, not significant in follow-up 12 weeks	no significant reduction	no significant reduction
<b>fatigue (MFIS)</b>	<b>significant reduction of fatigue</b> at end of sessions and in follow-up 2 weeks, gradual increase until follow-up 12 week	no significant reduction	no significant reduction

## Conclusions

The results show that HF-rTMS and iTBS significantly reduce spasticity, in contrast to placebo stimulation. Some evidence was found in favor of a longer-lasting effect of iTBS. A reduction in pain and fatigue could be achieved by HF-rTMS but did not result from iTBS.



# Abbreviations

<b>BS</b>	Brunnstrom stages
<b>CCAT</b>	Concise Chinese Aphasia Test
<b>CGI-S</b>	Clinical Global Impression–Severity scale
<b>DASS</b>	Depressions–Angst–Stress–Skala
<b>DLPFC</b>	dorsolateral prefrontal cortex
<b>FIM</b>	Functional Independence Measure
<b>FMA</b>	Fugl–Meyer Assessment
<b>GAD</b>	general anxiety disorder
<b>HAM-A</b>	Hamilton Anxiety Rating Scale
<b>HAM-D</b>	Hamilton Rating Scale for Depression
<b>HF-rTMS</b>	high-frequency rTMS
<b>iTBS</b>	intermittent theta-burst stimulation
<b>LANSS</b>	Leeds Assessment of Neuropathic Symptoms and Signs
<b>MADR</b>	Montgomery–Asberg Depression Rating Scale
<b>MAS</b>	Modified Ashworth Scale
<b>MBI</b>	Modified Barthel Index
<b>MFIS</b>	Modified Fatigue Impact Scale
<b>MSO</b>	Maximum Stimulant Output
<b>NAS</b>	Numerical Analog Scale
<b>NIHSS</b>	National Institutes of Health Stroke Scale
<b>OCD</b>	obsessive compulsive disorder
<b>SESS</b>	Subjective Evaluating Spasticity Scale
<b>SF-36</b>	Short Form Health Survey Questionnaire
<b>TPC</b>	Temporopolar Cortex
<b>U-FMA</b>	Fugl–Meyer Assessment for upper limb
<b>VAS</b>	Visual Analogue Scale
<b>VDS</b>	Verbal Descriptor Scale
<b>WMFT</b>	Wolf Motor Function Test
<b>Y-BOCS</b>	Yale–Brown obsessive compulsive scale

# Learn how to apply TMS

neurocare academy

Professionals who seek training through neurocare academy join a global network of professionals delivering best-practice neuromodulation. With online courses and practical workshops held each year throughout the world, the neurocare academy offers training in a range of applications including TMS, tDCS, advanced neurostimulation techniques, neurofeedback and applications of EEG.

Our new online learning platform means more practitioners from all over the world can access insights and know-how from a range of leading practitioners and scientists in the field. Professionals can register for instant access to learning modules and resources with the flexibility to complete a course at their own pace. This can then be followed by a practical workshop at one of our training centres worldwide.

We are a preferred training partner for private practices and hospitals worldwide and also offer custom in-house training and remote supervision for individuals and teams.

**For information,  
advice or to  
register:**

**neurocare group AG**

E-mail: [academy@neurocaregroup.com](mailto:academy@neurocaregroup.com)

phone: +49 (3677) 68 979-0

[www.neurocaregroup.com](http://www.neurocaregroup.com)

## Contact and information:

**neurocare group AG**

info@neurocaregroup.com

phone: +49 (3677) 68 979-0

[www.neurocaregroup.com](http://www.neurocaregroup.com)