

Accessories for the DC-**STIMULATOR** range

For our DC-**STIMULATOR** devices we offer a wide range of accessories, such as electrodes of various dimensions and shapes, cables and pastes. You can order these products individually or in a set. Please contact us. We will be happy to advise you which electrodes are best suited for your project.



Please send your order by e-mail to: sales@neurocaregroup.com

Content


Accessories for the DC- STIMULATOR range	1
Equipment starter set for DC- STIMULATOR / DC- STIMULATOR MOBILE	2
Equipment starter sets for DC- STIMULATOR MC	3
Rubber electrodes standard sizes (rectangular)	4
Sponge pads for rubber electrodes	5
Value packs of electrodes and sponges	6
Rubber electrodes suitable for MR	7
Rubber electrodes - special dimensions and shapes	8
Caps	10
Other materials (rubber straps, cables, pastes ...)	11
Contact for orders	12

Equipment starter set for DC-**STIMULATOR** / DC-**STIMULATOR** MOBILE

The equipment starter sets for the DC-**STIMULATOR** and the DC-**STIMULATOR** MOBILE consist of:

- 1 pair of rubber electrodes in different standard sizes (rectangular)
- 1 pair of sponge pads
- 1 pair of connection cables
- 1 combination of elastic rubber straps to attach the electrode to the head

All components can be ordered individually or in value packs (please see the following categories).



equipment starter set: 9 cm²
order # 305018



equipment starter set: 25 cm²
order # 305010



equipment starter set: 35 cm²
order # 305012



equipment starter set: 50 cm²
order # 305014



equipment starter set: 100 cm²
order # 305016

Equipment starter sets for DC-**STIMULATOR** MC

Each equipment starter set for the DC-**STIMULATOR** MC4 consists of:

- 4 pairs of rubber electrodes in different standard sizes (rectangular)
- 4 pairs of sponge pads
- 4 pairs of connection cables
- 4 combinations of elastic rubber straps to attach the electrode to the head

When using the DC-**STIMULATOR** MC8 or MC16 you require 2 or 4 starter sets.

All components can be ordered individually or as value pack (see following categories).



pack of electrodes MC4
3 x 3 cm (9 cm²)
order # 305026



pack of electrodes MC4
5 x 5 cm (25 cm²)
order # 305020



pack of electrodes MC4
5 x 7 cm (35 cm²)
order # 305021

Rubber electrodes standard sizes (rectangular)



3 x 3 cm, 9 cm²
1 pair
305054



5 x 5 cm, 25 cm²
1 pair
305050



5 x 7 cm, 35 cm²
1 pair
305051



5 x 10 cm, 50cm²
1 pair
305052



10x10cm, 100cm²
1 pair
305053

Sponge pads for rubber electrodes

The indicated dimensions are approximate values.



3 x 3 cm, 9 cm², 1 pair
305064



5 x 5 cm, 25cm², 1 pair
305060



5 x 7 cm, 35cm², 1 pair
305061



5 x 10 cm, 50 cm², 1 pair
305062



10x10 cm, 100 cm², 1 pair
305063



3 x 3 cm, 9 cm², 10 pairs
305064-11



5 x 5 cm, 25 cm², 10 pairs
305065



5 x 7 cm, 35 cm², 10 pairs
305066




5 x 10 cm, 50 cm², 10 pairs
305067




10x10 cm, 100 cm², 10 pairs
305068

Value packs of electrodes and sponges




value pack 9 cm²

order # 305112 (pack of 5)
order # 305122 (pack of 10)



value pack 25 cm²

order # 305113 (pack of 5)
order # 305123 (pack of 10)



value pack 35 cm²

order # 305114 (pack of 5)
order # 305124 (pack of 10)

Rubber electrodes suitable for MR

The indicated dimensions are approximate values.



suitable for MR
3 x 3 cm, 9 cm²
1 pair
order # 305059



suitable for MR
5 x 5 cm, 25 cm²
1 pair
order # 305055



suitable for MR
5 x 7 cm, 35 cm²
1 pair
order # 305056



suitable for MR
5 x 10 cm, 50 cm²
1 pair
order # 305057

Rubber electrodes - special dimensions and shapes

The indicated dimensions are approximate values.



Ø: 10 mm
area: 0,8cm²
without hole
2 mm thick
1 pair

order # 305090-01



Ø: 20 mm
area: 3 cm²
without hole
1 mm thick
1 pair

order # 305090-02



Ø: 25 mm
area: 5 cm²
without hole
2 mm thick
1 pair

order # 305090-11



Ø: 34 mm
area: 9 cm²
without hole
2mm thick
1 pair

order # 305090-03



Ø: 45 mm
area: 16 cm²
without hole
2 mm thick
1 pair

order # 305090-04



Ø: 75 mm
area: 44 cm²
without hole
2 mm thick
1 pair

order # 305090-05



Ø: 80 mm
area: 50 cm²
without hole
2 mm thick
1 pair

order # 305090-06



50 mm x 30 mm
area: 37 cm²
1 mm thick
1 pair

order # 305090-14



140mm x 70mm
area: 98 cm²
size: 2 mm
1 pair

order # 305090-10



outer Ø: 45 mm
 inner Ø: 15 mm
 area: 14 cm²
 with hole
 2mm thick
 1 pair

order # 305090-17



outer Ø: 45 mm
 inner Ø: 15 mm
 area: 14 cm²
 with hole
 2mm thick
 1 pair

order # 305090-17




outer Ø: 48 mm
 inner Ø: 24 mm
 area: 15 cm²
 with hole
 2 mm thick
 1 pair

order # 305090-12



outer Ø: 75 mm
 inner Ø: 20 mm
 area: 40 cm²
 with hole
 2 mm thick
 1 pair

order # 305090-07



outer Ø: 75 mm
 inner Ø: 30 mm
 area: 37 cm²
 with hole
 2 mm thick
 1 pair

order # 305090-15



outer Ø: 100mm
 inner Ø: 70mm
 area: 40 cm²
 with hole
 2 mm thick
 1 pair

order # 305090-16



outer Ø: 100mm
 inner Ø: 75 mm
 area: 40 cm²
 with hole
 2mm thick
 1 pair

order # 305090-08



outer Ø: 110mm
 inner Ø: 90mm
 area: 31 cm²
 with hole
 2 mm thick
 1pair

order # 305090-13

Caps



EEG electrode cap
no slots for adapters
chin strap

order # 905440



EEG electrode cap for ASIA
no slots for adapters
chin strap

order # 905441

Other materials (rubber straps, cables, pastes ...)



connection cable for rubber electrodes
approx. 150 cm, 1 pair (1-channel)
order # 305069



elastic rubber strap for rubber electrodes, 1 set
order # 305070



Extension rubber strap for rubber electrodes, 1 set
order # 305070-12



Velcro® strip for rubber electrodes, 1 set
order # 305072



extension of connection cables for rubber electrodes (5 m), 1 pair (1-channel)
order # 305075



extension of connection cables for rubber electrodes (10 m), 1 pair (1-channel)
order # 305076



Ten20 electrode paste
3 x 228 g cans
order # 905821



Ten20 electrode paste
3 x 114 g cans
order # 905822

Contact for enquiries and orders

neurocare group AG
Albert-Einstein-Straße 3
98693 Ilmenau

sales@neurocaregroup.com

+49 - 3677 – 68 97 90