

## Accessories and consumables for the THERA PRAX<sup>®</sup>

We offer a wide range of accessories for our THERA **PRAX**<sup>®</sup> devices, such as electrodes, pastes, caps or exercise materials. You can order these products individually, some are available in sets. Please contact us. We would be happy to advise you on which accessories are recommendable for your applications.

Please send your order by e-mail to: [sales@neurocaregroup.com](mailto:sales@neurocaregroup.com)

Please note that orders must have a minimum value of 75 EUR.



## Content

Accessories and consumables for the THERA <b>PRAX</b> <sup>®</sup> .....	1
Recommended accessories for clinical standard applications .....	2
Starter set for THERA <b>PRAX</b> <sup>®</sup> MOBILE and THERA <b>PRAX</b> <sup>®</sup> QEEG .....	4
Elektrodes and adapter .....	5
Consumables .....	6
Training materials.....	7
EEG caps.....	8
Biofeedback sensors .....	9
Contact for orders.....	10

## Recommended accessories for clinical standard applications

The following list contains accessories and consumables that are usually used in clinical neurofeedback trainings:



starter set for THERA **PRAX**® MOBILE  
(10 electrodes, Ten20, Nuprep, cotton buds)  
order # 905510



starter set for THERA **PRAX**® QEEG  
(10 electrodes, Ten20, Nuprep, cotton buds)  
order # 905515



1 sintered Ag/AgCl annular electrode  
(also suitable for ECG and EMG)  
order # 905105



10 sintered Ag/AgCl annular electrodes  
(also suitable for ECG and EMG)  
order # 905110



Nuprep abrasive paste (114g), pack of 3  
order # 905830



Ten20 electrode paste (114g can), pack of 3  
order # 905820



Ten20 electrode paste (228g can),  
pack of 3  
order # 905821



100 cotton buds, long, non-sterile  
order # 905920



sticking plaster, unsterile), 1 roll of 5 m  
order # 905940



electrode adapter for EEG caps (used with  
adhesive rings)  
order # 905200-10



500 adhesive rings for EOG conduction  
(20x8mm)  
order # 905910



"Neurofeedback" flyer for patient information  
order # 999910



transfer training cards for training at home  
order # 905550



DVD with video for transfer training at home  
order # 905570

## Starter set for THERA PRAX<sup>®</sup> MOBILE and THERA PRAX<sup>®</sup> QEEG



starter set of consumables for  
THERA PRAX<sup>®</sup> MOBILE

- 10 annular electrodes (also suitable for ECG and EMG),
- Ten20 electrode paste (1 x 114 g can)
- Nuprep abrasiv paste (1 x 114 g tube)
- 100 cotton buds, long

order # 905510



starter pack consumables for  
THERA PRAX<sup>®</sup> QEEG

- 23 annular electrodes (also suitable for ECG and EMG)
- Ten20 electrode paste (1 x 114 g can)
- Nuprep abrasiv paste (1 x 114 g tube)
- 100 cotton buds, long

order # 905515

## Elektrodes and adapter



sintered Ag/AgCl annular electrode  
(also suitable for ECG and EMG),  
1 unit

order # 905105



sintered Ag/AgCl annular electrodes  
(also suitable for ECG and EMG),  
10 units

order # 905110




electrode adapter for EEG caps (with  
905910 adhesive rings, see consumables)

order # 905200-10



snap cable for disposable EMG/ECG  
electrodes (150 cm), 10 pieces

order # 905120



ECG electrodes, snap, disposable, 50  
pieces

order # 905121

## Consumables



Nuprep abrasive paste (114g), pack of 3  
order # 905830



Ten20 electrode paste (114 g can),  
pack of 3  
order # 905820



Ten20 electrode paste (228g can),  
pack of 3  
order # 905821



500 adhesive rings for EOG conduction  
(20x8mm), (with 905200-10 electrode adapter)  
order # 905910



Sticking plaster, unsterile, 1 roll of 5 m  
order # 905940



cotton buds, long, unsterile, 100 pieces  
order # 905920



electrode gel ABRALYT HiCl (1,000 g)  
order # 905810



EEG gel OneStep (120 g)  
order # 905825

## Training materials



transfer training cards for training at home  
order # 905550



transfer training video, for training at home  
order # 905570

## EEG caps

All caps are available in the sizes 52, 54, 56, 58 and 60. Please indicate the size when ordering.



EEG electrode cap TP, 14+4 additional adapters, chin strap

order # 905240



EEG electrode cap TPQ-EEG, 23+4 additional adapters, chin strap

order # 905245



Starter set for EEG caps (ABRALYT HiCl, syringes, cotton buds, adaptors, sticking rings, Ten20, Nuprep)

order # 905800



## Biofeedback sensors

For biofeedback applications your device requires the 4-channel biofeedback extension. Please contact us. We will be happy to submit you a quotation:

[info@neurocaregroup.com](mailto:info@neurocaregroup.com)



pulse curve sensor  
reflection sensor

order # 105010



pulse curve sensor  
ear clip

order # 105011



pulse curve sensor  
finger clip

order # 105012



GSR sensor  
(galvanic skin response)

order # 105020



breath sensor

order # 105030



temperature sensor

order # 105040

## Contact for enquiries and orders

neurocare group AG  
Albert-Einstein-Straße 3  
98693 Ilmenau

[sales@neurocaregroup.com](mailto:sales@neurocaregroup.com)

+49 - 3677 - 68 97 90